

Program start

ema2wav can be started either

- ▶ using the **GUI**
 - ▶ via binary (.dmg)
 - ▶ via console
- ▶ using command-line (see documentation)
- ▶ by importing as a python module (custom script/notebook/**Google colab notebook**; see documentation)

Installation (Anaconda - script version)

- ▶ open Terminal/Console/Anaconda Prompt
- ▶ create conda environment:
`conda create -name ema_env`
- ▶ activate the environment:
`conda activate ema_env`
- ▶ install the dependencies:
`pip install -r`
`path/to/ema2wav/folder/src/requirements.txt`

Run ema2wav

- ▶ for Mac: run `ema2wav_app.app` (see GitHub documentation)
- ▶ script version (GUI, from terminal, `ema2wav src` directory assumed):
`python ema2wav_app.py`


User input

```
{  
  "ema_device_info": "AG50x",  
  "export_audio+ema": true,  
  "export_to_csv": false,  
  "export_raw_ema": false,  
  "output_directory": "/path/to/your/output/folder/",  
  "ema_input_directory": "/path/to/your/ema/folder/",  
  "ema_samplerate": 1250,  
  "ema_channels": 16,  
  "audio_input_directory": "/path/to/your/wav/folder/",  
  "audio_samplerate": 48000,  
  "audio_channels": 1,  
  "channel_allocation": {  
    "TD": 3,  
    "TT": 5  
  },  
  "parameters_of_interest": {  
    "0_TT": "y",  
    "1_TD": "y-vel",  
    "2_TT": "tvel",  
    "3_TD": "x"  
  },  
  "filter": null  
}
```

User input - parts of the input file

<code>{</code>		
<code>"ema_device_info": "AG50x",</code>	Device info (necessary)	
<code>"export_audio+ema": true,</code>		
<code>"export_to_csv": false,</code>	output options (necessary)	
<code>"export_raw_ema": false,</code>		
<code>"output_directory": "/path/to/your/output/folder/",</code>		
<code>"ema_input_directory": "/path/to/your/ema/folder/",</code>	folder paths (necessary)	
<code>"audio_input_directory": "/path/to/your/wav/folder/",</code>		
<code>"ema_samplerate": 1250,</code>		
<code>"ema_channels": 16,</code>	general information of the ema	
<code>"audio_samplerate": 48000,</code>	and audio files (optional)	
<code>"audio_channels": 1,</code>		
<code>"channel_allocation": {</code>		
<code> "TD": 3,</code>	channel information	
<code> "TT": 5</code>	(necessary)	
<code>},</code>		
<code>parameters_of_interest": {</code>		
<code> "0_TT": "y",</code>	parameters of interest,	
<code> "1_TD": "y-vel",</code>	consist of X_ + channel name	
<code> "2_TT": "tvel",</code>	& the parameter to extract	
<code> "3_TD": "x"</code>	(necessary)	
<code>},</code>		
<code>"filter": null</code>	filter information (necessary)	
<code>}</code>		

Start of the conversion process (GUI)




Start of the conversion process (manual)

- ▶ from terminal (ema2wav src directory assumed)
`python convert.py config.json`
- ▶ as Python module:

```
import ema2wav_core as ec  
config_file =  
"/path/to/your/config_file.json"  
ec.ema2wav_conversion(config_file)
```

Conversion process



Conversion process - reshaping of the data

sample 1				sample 2				sample n
channel 1				channel 1				...
x	y	z	...	x	y	z




sample n							
	x	y	z	phi	theta	rms	empty
channel 1
channel 2
channel n

sample 2							
	x	y	z	phi	theta	rms	empty
channel 1	value	value	value	value	value	value	value
channel 2	value	value	value	value	value	value	value
channel n


sample 1							
	x	y	z	phi	theta	rms	empty
channel 1	value	value	value	value	value	value	value
channel 2	value	value	value	value	value	value	value
channel n

Conversion process - Interpolation



(https://en.wikipedia.org/wiki/Spline_interpolation/media/File:Cubic_splines_three_points.svg)

Conversion process - GUI approach



Conversion process - GUI approach (2)

- ▶ open the folders containing .pos and .wav files
- ▶ channel allocation: enter name channel number
- ▶ parameters of interest: enter channel name parameters
- ▶ select filter (if necessary)
- ▶ open output folder
- ▶ start the conversion

Conversion process - GUI approach (3)

- ▶ config file as documentation
- ▶ can be used for replicating the conversion process (load CONFIG)

Conversion process - Google Colab approach

- ▶ **GOOGLE COLAB**
- ▶ execute Jupyter Notebooks online
- ▶ free easy-to-use
- ▶ see example notebook

Conversion process - Google Colab (1)

- ▶ open the example notebook in Google Colab

[illegible]

Conversion process - Google Colab (2)

- ▶ open the files folder

[illegible]


Conversion process - Google Colab (3)

- ▶ run the first code snippet

[illegible]

Conversion process - Google Colab (4)

- ▶ upload your .pos and .wav files




Conversion process - Google Colab (5)

- ▶ upload the GC config file

[illegible]

Conversion process - Google Colab (6)

- ▶ run the other code snippets and download converted files



```
1. upload your wav and pos files into the input folder
2. upload the config file

import ena2wav

[2]: import requests

      https://raw.githubusercontent.com/gh00bz/research_wav2vec2_resnet50_pytorch/main/ena2wav.py
      requests.get(url)
      wvc = wav2vec2_pytorch.Wav2Vec2Wrapper.from_instances([wvc])
      wvc.save_pretrained('wvc')

[3]: !pip install Outage

      Looking in indexes: https://pypi.org/simple, https://www.pycharm.org/pypi, https://mirrors.aliyun.com/pypi/simple/
      Collecting Outage
      Downloading Outage-1.45.1-py3-none-any.whl (228 kB)
      Installing collected packages: Outage
      Successfully installed Outage-1.45.1

import ena2wav

[4]: import Outage as kb

      Start conversion process


      kb.wav2vec2_resnet50_pytorch('wvc')
```

Annotation & Measurements in Praat

Praat tweaks

- ▶ Display default settings in Praat are not suitable for annotating EMA data
- ▶ Options for the best annotation experience:
 - ▶ disable the spectrogram
 - ▶ Change sound scaling in the editor window:
(Sound > Sound scaling... > select 'by window and channel')
 - ▶ Mute channels (if necessary):
(Sound > Mute channels... > (ranges) > enter 2:X)


Landmark annotations



Gestural landmarks and intragestural intervals for a typical gestural complex (Tilsen, 2014)

Annotation & Measurements in Praat

Annotations



- ▶ Annotations of gestural landmarks using Point Tiers in Praat
- ▶ retrieve timing information as usual
- ▶ measure amplitude information by calling 'Get value at time' applied to the Sound object
- ▶ see **measurement_example.praat**

Remarks

- ▶ never use the entire file for frequency/intensity/spectral measurements!
Instead: extract the the audio track and use that for acoustic measurements
- ▶ make use of a good documentation of your files